

# Safety card

## Pack size

### Designation

SAFETY CARD AB

SAFETY PACK

1 Disposable lancet + 1 securitube + 1 vial of isotonic saline solution + 4 strirrer-pipettes

SAFETY TEST ABO

1 SAFETY CARD AB + 1 SAFETY PACK

### Reference

79074

79055


79069

### quantity


100 cards

10 kits


100 kits



## Rules of compatibility




## Exemples of results




Reactions result in agreement. According to the ABO compatibility rules, transfusion is possible

**Compatible**



Reactions result in disagreement. According to the ABO compatibility rules, transfusion is possible

**Compatible**



Reactions result in disagreement. According to the ABO compatibility rules :

**transfusion is not possible**

**Incompatible**

## Customization

Bedside tests are available in your language. Feel free to contact us for further information.

