


Accessory Sets

OriGen supplies a variety of spike adapters for cell transfer, recovery, and administration. All are packed individually in a peel pouch and are sold sterile. Case quantities vary, see the description below. All materials are USP Class IV devices.



Spike Adapters

Product Code	Description	Case Quantity	Stock or Custom
CBS	Needle free adapter port on a spike adapter, swabbable with alcohol, 100% DMSO resistant	50	Stock
IBS	Standard needle injection port on a spike adapter, swabbable with alcohol, 100% DMSO resistant	50	Stock
S-38B12	Spike Adapter x 3/8 Barb	25	Custom

Empty Syringes

Product Code	Description	Case Quantity	Stock or Custom
RF-T15	6ml syringe with female luer and 15cm SCD tubing	25	Stock
RV-T15	6 ml syringe, Needle free adapter with 15cm SCD tubing	25	Custom
RV-T50	50ml syringe with needle free adapter and 15cm SCD tubing	20	Custom

Other Accessories

NFV	Needle-Free Adapter male luer	10	Stock
CT	CryoTape - Stable in LN temperatures	10	Stock

Accessory Sets



Tube Sets

Product Code	Description	Case Quantity	Stock or Custom
4S-4M60	Manifold Set: Four spike adapters to four male luers with 60cm of PVC tubing and pinch clamps	12	Stock
S-F10	Spike adapter with female luer, 10cm PVC/EVA tube and pinch clamp, 100% DMSO resistant	50	Stock
S-M10	Spike adapter with male luer, 10cm PVC/EVA tube and pinch clamp, 100% DMSO resistant	50	Stock
S-F50	Spike adapter with female luer, 50cm PVC tube and pinch clamp	50	Custom
S-M50	Spike adapter with male luer, 50cm PVC tube and pinch clamp	50	Custom
S-S50	Spike adapter to spike adapter, 50cm PVC tube and pinch clamp	50	Custom
2P-S16	One spike adapter to two spike ports, 16cm PVC tubing	20	Custom
F-M12	Female luer to male luer, 12cm of PVC tubing	50	Custom
SU-M50	Spike adapter to male luer adapter, 50cm PVC tubing, roller clamp	50	Custom

All spike adapters, tube sets are resistant to 80% DMSO unless otherwise specified.

Other custom configurations available. Minimum quantities may apply.

OriGen
BIOMEDICAL

Americas & Asia

7000 Burluson Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Fax: Toll Free (888) 812 8411 or
+1 (512) 617 1503

Europe & Middle East

Landskronavägen 25A
252 32 Helsingborg, Sweden
Tel. + 46 46 772 433
Fax. +46 44 785 1099